



Treatment Results



Courtesy of Dr. Adina Margaritescu,
Bucharest, Romania



Courtesy of Xavier Martens, MD, Euroclinix,
Hasselt, Belgium




Courtesy of Yuichi Sal, MD,
Japan



Courtesy of Dr. Rasya Dixit,
Bangalore, India

Novoxel®



Exclusive Distributor:

REDO
MARKETING INDONESIA
Your Aesthetic Business Solutions Provider

Jalur Sutera Timur, Ruko Renata Blok 1A/9
Alam Sutera, Tangerang 15144
Tel: 021-2966 1469 Fax: 021-3043 9750
Website: www.redo.co.id
Email: enquiry@redo.co.id
All Right Reserved by PT. Redo Marketing Indonesia

What Is Tixel®?

- A novel non-ablative & ablative fractional skin rejuvenation system powered by Novoxel's patented titanium tip technology.
- Tixel treats delicate facial skin, including periorbital, eyelids, the neck and decollete safely and quickly with low pain and 0-2 days downtime.
- Tixel employs direct heat and is radiation-free.



MULTIPLE TREATMENTS
Skin rejuvenation, crow's feet, open pores, acne, stretch marks...



MULTIPLE AREAS
Face, neck, decollete hands, arms & body



LOW PAIN
Analgesic creams are not required



STERILE TITANIUM TIP
The Tip biocompatible and sterile during treatment



SMALL & PORTABLE
Small, light, easy to carry



RADIATION FREE
Tixel employs direct heat



Indications

- Wrinkles and fine lines
- Skin discoloration
- Sun damaged skin
- Skin tone and texture
- Tightening of the upper eyelids
- Age spots
- Stretch marks / Striae

Thermal energy is transferred to the skin by a hot titanium Tip, which contacts the skin briefly. There is no bleeding.

The extent of energy transfer is determined by the Tip's exposure time with the skin.

Heat transfer is by contact. The Tip does not perforate the tissue.



Forehead
"worry" lines

Glabellar
"frown" lines

Eyelid lines

Acne scars

Lip lines

Creased chin

Crow's feet
(Periorbital lines)

Bunny lines

Nasal labial folds
(Smile lines)

Neck folds

